


Measurements of ^{223}Ra

Jeffrey T. Cessna


*Physics Laboratory, National Institute of Standards and Technology,
Gaithersburg, MD 20899-20899, USA*

Motivation

- Nuclear Medicine Standardization
- ^{223}Ra used in the treatment of skeletal metastases and for therapy of associated bone pain
- Naturally bone-seeking
- In clinical trials for treatment of skeletal metastases from prostate cancer and breast cancer
- Food and Drug Administration required a NIST standard before initiation of clinical trials in the US



^{99m}Tc -MDP



^{223}Ra

Anterior

Posterior

Figure 4. Scintigraphic images demonstrating accumulation of ^{223}Ra in skeletal lesions in accordance with ^{99m}Tc -MDP uptake. Radium image taken 24 hrs post injection, hence radium excreted and in transit in the large bowel is evident in the lower left image. (Image quality is dependent on dose administered: 750 MBq ^{99m}Tc -MDP vs. 12 MBq ^{223}Ra)

ECCO12 Poster No. 543
Copenhagen, 2003

Experiments

- E1 – Preliminary measurements
- E2 – Adsorption studies
- E3 – Check for other possible losses, preliminary $2\pi\alpha$ proportional counting (PC)
- E4 – LS counting including PTB composition, $2\pi\alpha$ PC
- E5 – Final LS counting, Final $2\pi\alpha$ PC
- E6 – Dose calibrator settings for vial

E3 - $2\pi\alpha$ PC

- Try $2\pi\alpha$ gas-flow proportional counting as a second confirmatory measurement, observes only alphas
- 1.5 cm stainless steel disks, 1-2 drops active
- Learn we should cover the sources with collodion film
- Lynne cleans the detector

results agree with LS measurements

Experiment E4

Stock solution received from Algeta, ASA

E4: 20 MBq in 10 mL

Solution consisted of 1.5 mL ^{223}Ra solution in carrier*, diluted to 10 mL with 8.5 mL of 0.9 mol·L⁻¹ NaCl solution containing 28 mmol·L⁻¹ Na₃C₆H₅O₇.

Counting in NIST "4π"γ ionization chamber (+ activity calibrators)

Transfer nominally 5 mL to each ampoule

E4A1

E4A2

Dilute A1 (DF = 3.944) with carrier*

Bring remainder volume to 5 mL with carrier*

Dilute A2 (DF = 3.982) with 6 mol·L⁻¹ HNO₃

E4A1D1

E4A1R1

E4A2D1

E4A2R1

Counting in NIST "4π"γ ionization chamber, activity calibrators

Dilute A1D1 (DF = 50.263) with carrier*

E4A1D2

Liquid scintillation counting (LSC) with ^3H -standard efficiency tracing

Measurement on HPGe systems

Counting in NIST 2π α proportional counter

*Carrier consisted of nominally 5 mmol·L⁻¹ Na₃C₆H₅O₇, 2 mmol·L⁻¹ Sr²⁺, and 0.5 mmol·L⁻¹ $^{45}\text{Ca}^{2+}$.

Experiment E5

Counting in NIST “4 π ” γ
ionization chamber



Stock solution received from Algeta, ASA

E5: 10 MBq in 5 mL

Solution consisted of 0.75 mL ^{223}Ra solution in carrier*,
diluted to 10 mL with 4.25 mL of 0.9 mol·L⁻¹ NaCl solution
containing 28 mmol·L⁻¹ Na₃C₆H₅O₇.



*Transfer
contents and
bring to 5 mL
with carrier**

E5A1

*Dilute A1 (DF = 4.864) with
carrier**




E5A1D1

*Dilute A1D1 (DF = 46.042) with
carrier**



E4A1D2

Counting in NIST 2 $\pi\alpha$
proportional counter



Counting in NIST “4 π ” γ
ionization chamber, activity
calibrators

Liquid scintillation counting (LSC) with
 ^3H -standard efficiency tracing

Measurement on
HPGe systems



*Carrier consisted of nominally 5 mmol·L⁻¹ Na₃C₆H₅O₇, 2 mmol·L⁻¹ Sr²⁺, and 0.5 mmol·L⁻¹ Ca²⁺.

E4 & E5– $2\pi\alpha$ PC

- Repeat $2\pi\alpha$ PC
- Source preparation
 - 1.5 cm stainless steel disks
 - 1-2 drops active material
 - collodion film covers to 2 to 8 $\mu\text{g}\cdot\text{cm}^{-2}$
- Weigh a group of covers to determine average thickness
- Perform extrapolation to no covering
 - Using determined thickness
 - Using number of covers
 - Agreement between extrapolations

E6 – Dial Settings

- 5 ml NIST ampoule, 5 ml solution
- 20 ml FIOLEX (MGlaser AG, Mönnerstadt, Germany) dose vial
 - 3 x 0.5 ml solution
 - 3 x 2 ml solution
 - 3 x 6 ml solution
- Radionuclide Calibrators
 - Capintec CRC-12, 15R, 35R
 - Biodex Atomlab 100
 - NPL 671

E6 – Preliminary Results

	Capintec CRC-12		Capintec CRC-15R		Capintec CRC-35R		Atomlab 100		NPL Chamber		
	Dial Setting	U _c	Dial Setting	U _c	Dial Setting	U _c	Dial Setting	U _c	Calibration Factor		U _c
5 mL ampoule	264	4	265	4	266	4	16.9	0.2	3.49899	pA/MBq	0.04506
20 mL Dose vial	261	4	263	4	266	4	16.9	0.2	3.50187	pA/MBq	0.04504